

Software License Agreement

This is a legal agreement between you (hereinafter referred to as LICENSEE) and Bruker Optics GmbH & Co. KG, Rudolf-Plank-Straße 27, 76275 Ettlingen, registered in the commercial register of the Local Court of Mannheim under HRA 710333 (hereinafter referred to as LICENSOR). Carefully read all terms and conditions of this agreement before opening the sealed disc package. Opening the sealed disc package indicates your acceptance of the following terms and conditions of this agreement. If you do not agree to the terms and conditions of this agreement, return the unopened disc package to the LICENSOR.

1. Subject of agreement

The program OPUS and software extensions of OPUS (hereinafter referred to as SOFTWARE) and copyrighted libraries and other copyrighted digital data stored on the delivered data medium (hereinafter referred to as DATA) and the accompanying documentation are subject of this license agreement. The LICENSOR points out that, according to the state of the art, it is not possible to create SOFTWARE and/or DATA in a way, that they operate error free in any application and combination. Therefore, only the SOFTWARE and/or DATA usable in terms of the user manual and product description are subject of this agreement.

2. Grant of license and restrictions

The LICENSEE may install, run and/or use the SOFTWARE and/or DATA stored on the delivered data medium on one (1) computer. If the LICENSEE wants to use the SOFTWARE and/or DATA on more than one (1) computer at the same time the LICENSEE needs to purchase additional licenses.

The LICENSEE may store or install copies of the SOFTWARE and/or DATA on networks server only on the condition that this copy is used to install or run the SOFTWARE and/or use the DATA via network on one (1) computer. Sharing or using one license for several computers at the same time is strictly forbidden. If LICENSEE is using OPUS in combination with the OPUS/DLL package OPUS can be installed on a second computer to allow configuration and servicing of the OPUS/DLL installation.

3. Transfer restrictions

The LICENSEE may not transfer, loan, rent, lease or sell the SOFTWARE and/or DATA or parts thereof either free of charge or against payment to any third party. If a transfer of the SOFTWARE and/or DATA license is unavoidable (e.g. by selling the spectrometer) this transfer can only be done by the LICENSOR. The LICENSEE is obligated to inform the LICENSOR in time and in writing about the intended transfer.

4. Other restrictions


The LICENSEE may not modify, adapt, translate, reverse engineer, decompile or disassemble the SOFTWARE in order to obtain the source code. The LICENSEE may not use the DATA in a manner other than stated in the product documentation or extract information from DATA in whatever form. Printouts or reproductions of the DATA or parts thereof must be provided with a copyright note.

5. Copyright and proprietary rights

SOFTWARE, DATA, the accompanying documentation and each copy of SOFTWARE, DATA and accompanying documentation are proprietary of the LICENSOR and are protected by copyright laws, international treaty provisions and all other applicable national laws. Therefore, the LICENSEE must treat the SOFTWARE, DATA and accompanying documentation like any other copyrighted material except that the LICENSEE may either (a) make one (1) copy of the SOFTWARE and/or DATA solely for backup or archival purposes or (b) transfer the SOFTWARE and/or DATA to a single hard disk provided the LICENSEE keeps the original solely for backup or archival purposes. The LICENSEE may not copy the written materials accompanying the SOFTWARE and/or DATA. The LICENSOR points out that the LICENSEE is liable for any damages arising from any copyright infringement.

6. Warranty and limitation of liability

The LICENSOR warrants that the SOFTWARE and/or DATA will perform substantially in accordance with the accompanying product manuals for a period of 24 months from the date of receipt. Furthermore, the LICENSOR warrants that the data medium on which the SOFTWARE and/or DATA are recorded shall be free from defects in materials and workmanship under normal use and service for the same period. If the data medium, on which SOFTWARE and/or DATA are recorded, is defective, the LICENSEE must return the defective data medium to the LICENSOR and the LICENSOR will replace it at no charge.



For the reasons named under provision 1 the LICENSOR does not warrant that the SOFTWARE, DATA and accompanying documentation are error free. The LICENSOR disclaims all warranties, either express or implied, including but not limited to warranties of fitness for a particular purpose and the compatibility with other software programs or hardware components. The LICENSEE assumes full responsibility for the selection of the software and the consequences resulting from the use of the software. The LICENSEE assumes also the entire risk as to the intended or achieved results. The LICENSOR is only liable in cases of intention and gross negligence.

The LICENSOR shall not be liable for any special, incidental, indirect or consequential damage of any kind or other damages whatsoever (including, but not limited to damages for loss of business profits, business interruptions, loss of business information or other pecuniary loss) arising out of or in connection with the use or non-use of the SOFTWARE and/or DATA, even if the LICENSOR has been notified of the possibility of such damages. In any case the LICENSOR'S entire liability shall not exceed the actual amount of money paid by the LICENSEE for SOFTWARE and DATA. Some states and jurisdictions do not allow the exclusion or limitation of liability for consequential or incidental damages, so the above limitation or exclusion may not apply in any case.

The SOFTWARE is pre-configured such that the LICENSEE can use it in compliance with data protection. The LICENSEE is, however, responsible for the effective implementation of the software use in compliance with data protection. In this connection, the LICENSOR does not warrant that the LICENSEE has implemented all data protection requirements of the European General Data Protection Regulation (GDPR) when using the software. Any claims made against the LICENSOR which are based on the violation of data protection obligations on behalf of the LICENSEE are excluded.

7. Term and termination

This license is effective until terminated. The LICENSEE'S rights under this license will terminate automatically without notice from the LICENSOR if the LICENSEE fails to comply with any term or condition of this license. In this case the LICENSEE is obligated to return immediately the original and all copies of SOFTWARE, DATA and accompanying material to the LICENSOR.



8. Miscellaneous provisions

Place of jurisdiction for any legal dispute arising out of or relating to this agreement is Karlsruhe. If any provision of this agreement is deemed invalid or incomplete, partly or entirely, the remaining provisions will remain in full force and effect.

Status: December 2021



OPUS

Introduction

© 2021 Bruker Optics GmbH & Co. KG, Rudolf-Plank-Straße 27, D-76275 Ettlingen
www.bruker.com

All rights reserved. No part of this user manual may be reproduced or transmitted in any form or by any means including printing, photocopying, microfilm, electronic systems etc. without our prior written permission. Brand names, registered trademarks etc. used in this user manual, even if not explicitly marked as such, are not to be considered unprotected by trademarks law. They are the property of their respective owner.

The description given in this user manual is based on the technical specifications and the technical design valid at the time of publication. Technical specifications and technical design may be subject to change.

This manual is the original documentation for the OPUS spectroscopy software.

Table of Contents

1	Introduction	5
1.1	About this manual	5
1.2	Terms.....	5
1.3	Gender-neutral form.....	5
1.4	Questions and concerns.....	6
2	Installation Instructions.....	7
2.1	System requirements for software	7
2.2	Access rights.....	9
2.3	Miscellaneous	9
2.4	Installing OPUS using Windows 10	9
2.5	Modifying, repairing or removing single OPUS features	10
3	Installing OPUS	11
3.1	Installation procedure	11
3.2	OPUS Login.....	12
3.2.1	Changing password in the validation mode	13
3.3	OPUS user interface.....	14
3.4	Access control.....	17
3.4.1	User record	17
3.4.2	OPUS workspace	18
3.4.3	Data recording	18

Table of Contents

4	Working with OPUS	19
4.1	Measuring a spectrum.....	19
4.1.1	Connecting to spectrometer.....	19
4.1.2	Setting up measurement parameters	21
4.1.3	Calibrating interferogram peak	24
4.1.4	Storing experiment file	25
4.1.5	Measuring background spectrum	27
4.1.6	Measuring sample spectrum.....	27
4.2	Baseline correction	30
4.2.1	Performing baseline correction	31
4.3	History of spectrum files.....	33
4.4	Viewing global audit trail	35
4.4.1	Loading HISTORY data block from global audit trail	36
4.5	Saving spectrum	36
4.5.1	Saving and overwriting the original file	36
4.5.2	Saving without overwriting the original file	37
4.6	Load spectrum file.....	38
4.6.1	Description of Load Spectrum dialog.....	39
4.7	Printing spectrum file.....	41
4.7.1	Using the Quick Print command	41
4.7.2	Using the Print Spectra command	41
5	Analysis Report.....	43
5.1	Generating report	43
	Index.....	45

1 Introduction

1.1 About this manual

This manual is a quick reference guide about the OPUS spectroscopy software. For important information the manual uses the following format:

i The *i* provides important information given to the user to facilitate system operation, and make the best possible use of the software.

1.2 Terms

This manual uses the term *OPUS* when the software is described.

1.3 Gender-neutral form

The user manual uses the male form in a neutral sense and does not differentiate between male, female and diverse users. We would kindly ask the users to have understanding for the male form which only serves to simplify the creation of this manual.

1.4 Questions and concerns

If you have questions or concerns about operating the software, contact Bruker at the numbers listed below:

- Service hotline: +49 (0) 72 43 504-2020
- Fax: +49 (0) 72 43 504-2100
- E-mail: service.bopt.de@bruker.com
service.bopt.us@bruker.com
- Internet: www.bruker.com/optic
- International service: www.bruker.com/about-us/offices/offices/bruker-optics

2 Installation Instructions

2.1 System requirements for software

PC data system:	<ul style="list-style-type: none">• serial interface^a• Ethernet, RJ45, 10Base T (10 or 10/100MHz)• USB^b• Mouse, keyboard• PCI-x1 interface^c
Processor:	64 bit processor
RAM:	4 GB
Hard disk:	250 GB, IDE
Graphics card:	XVGA, >128MB RAM, 32 color-bit depth, open GL and support of Direct-X, version 9 or higher (for OPUS/VIDEO, OPUS/3D)
Monitor:	Minimum resolution: 1024 x 768 Recommended resolution: 1600 x 1200
Operating system:	Minimum (OPUS 8 and higher): Windows 10 (Home Premium or better 64 bit)

Table 2.1: System requirements for software

Ethernet:

Windows 10 with OPUS 8 or higher:

- Internet explorer 11

General (for any OPUS version):

- TCP/IP must be configured
- Address range of PC network card must include spectrometer IP
- Port 80 to spectrometer via network card must be available
- Proxy is not allowed to interfere with the communication to the optics
- No USB port restrictions

Table 2.1: System requirements for software

- a. Only required when connecting external devices, e.g. temperature controller, xy stage controller etc.
- b. USB 2.0 is required to load OPUS software, and operate the digital video cameras built-in the Ramanscope III, SENTERRA (2x). USB 3.0 is required to operate the LUMOS IR microscope.
- c. Some types of accessories (e.g. A673-X, A750/x, S522-20, S522-200, Frame-Grabber etc.) require a data system which is able to accommodate one or more cards with full length.

2.2 Access rights

- Windows Administrator rights required to install or register OPUS
- Complete write access required to the target path
- Complete access required to USB ports, no restrictions on access rights by services such as *DeviceLock Service*

2.3 Miscellaneous

Configure Windows so that it is possible to see the file extensions of registered file types.

1. Open the Windows Explorer.
2. On the *Tools* menu, click the *Folder Options* command.
3. Click the *View* tab.
4. Delete the checkmark in front of the *Hide extensions for known file types* option.
5. Click *OK* to confirm the settings.

2.4 Installing OPUS using Windows 10

When using Windows 10 it is not sufficient to be logged into the operating system as a user with Administrator rights. You need to perform the installation explicitly as administrator (chapter 3.1).

2.5 Modifying, repairing or removing single OPUS features

i For technical reasons not all features are deleted when de-installing the OPUS program. If you want to remove OPUS completely, you have to delete the OPUS directory manually after de-installation. Before you do this, store all important data (e.g. workspaces, printing templates, macros etc.) in a different directory.

1. Click the Windows Start menu.
2. Select the *Control Panel* command.
3. On the view that opens, select the *Programs and Features* command.
4. Right click the current OPUS version installed.
5. From the pop-up menu displayed, select the *Change* command.
6. On the view that opens, check the particular option button depending whether you want to modify, repair or remove OPUS features.
7. Click the *Next* button.
8. Use the InstallShield Wizard to make the settings.
9. Click the *Finish* button.

3 Installing OPUS

i By default, the OPUS spectroscopy software is supplied on a USB stick.

3.1 Installation procedure

i Installing OPUS requires Administrator user rights.

1. Insert the OPUS stick supplied into the USB port of your PC.
2. If the installation starts automatically after inserting the OPUS stick, interrupt the installation procedure.
3. In the *OPUS_<version>* installation directory, right click the *setup<nn>.exe* file.
4. From the pop-up menu that opens, select the *Run as Administrator* option.
5. Confirm the message that pops up to be able to continue.
6. Follow the on-screen installation instructions and enter the OPUS registration data.¹
 - *A print-out of the registration data is annexed to the OPUS stick. The registration has to be performed only once and is valid for all users.*

1. If you want to install an OPUS test version (valid for 30 days), you do not need to enter any registration data. In this case, set a checkmark for the test option.

7. If OPUS is to run in the 21 CFR Part 11 operating mode (validation mode), activate the particular option.¹
 - *If OPUS 8 or higher is to run in the validation mode, it is recommended to work with a database. The connecting option to a database (OPUS database or user-specific database) can be activated during installation. If required, the setting made can only be undone later by a user with Windows administrator rights.*
8. Select the language to be used for OPUS.
9. Restart the PC at the end of the installation procedure.
 - *If you are prompted to restart your computer when you install different program components, do NOT restart before the OPUS setup has been completed.*

3.2 OPUS Login

1. When the installation has been finished, click the OPUS icon on the Windows Start menu.
2. On the dialog that opens, select the user ID² from the drop-down list.
3. Enter *OPUS* as password.
 - *The password is case sensitive. Once you have assigned yourself a user account in OPUS, it is possible to determine your own user ID and password. When working in validation mode, you will be requested to change the password immediately after the first login. Details are described in chapter 3.2.1.*
4. Select a workspace.

-
1. The option can only be activated if the OPUS registration contains the OPUS/VALIDATION software package. Activating the option cannot be undone and has comprehensive effects when working with OPUS. 21 CFR part 11 is a quality regulation legally approved by the United States Food and Drug Administration (FDA), which also regulates the requirements on data integrity. If OPUS 8 or higher is to run in the validation mode, it is recommended to work with a database. The option to be connected to a database (OPUS database or user-specific database) can only be activated during installation.
 2. By default, OPUS provides the following user IDs: Admin, LabManager and Operator. Detailed information on the user ID are described in the OPUS Reference Manual.

- *You can always select those workspaces from the drop-down list which have been defined for your user ID. When you use OPUS for the very first time, we recommend to accept the standard 'default.ows' workspace.*
- 5. Click the *Login* button.
- 6. On the dialog that opens, click the *OK* button.
 - *The OPUS user interface opens.*

3.2.1 Changing password in the validation mode

1. After the very first OPUS login you are requested to change the password. Click *OK* to confirm the request.
2. In the dialog that opens, enter the new password and verify it.
 - *According to the 21 CFR Part 11 operating mode (validation mode) the password has to comply with special complexity requirements, therefore has to fulfill the following minimum conditions:*
 - *Total length: 8*
 - *Number of numerals: 1*
 - *Number of letters: 1*
 - *Number of lower case letters: 1*
 - *Number of upper case letters: 1*
 - *Number of special characters: 1*
3. Click the *Change Password* button.
4. In the *Login* dialog, enter the new password and click the *Login* button.

3.3 OPUS user interface




Figure 3.1: OPUS user interface

	Element	Definition
1	Title bar	The title bar shows the user ID and the workspace assigned to the user.
2	Menu bar	The menu bar shows all menus available in OPUS. The individual menus allow access to all OPUS commands.

Table 3.2: Elements of OPUS user interface

	Element	Definition
3	Toolbar	The toolbar contains frequently-used OPUS commands, individually configured by means of icons.
4	Browser	The browser shows the spectrum files loaded, with the particular data blocks.
5	Spectrum window	The spectrum window shows the spectrum file loaded. By default, a crosshair is displayed. The current position of the crosshair (e.g. wavenumber and intensity) are displayed on the upper right edge. This allows to read out the x/y positions directly from the spectrum. Right clicking the spectrum window hides the crosshair.
6	Overview window	The overview window shows the complete frequency range of the spectrum file loaded.
7	OPUS view	The OPUS view you work with is indicated by a tab. Each tab has a different color to be able to distinguish between the OPUS views. If several OPUS views are open simultaneously, click the respective tab or use the buttons on the left and right to change between the different OPUS views.
8	Online help	The OPUS online help provides detailed information on all OPUS commands and can be opened by clicking the F1 key.

Table 3.2: Elements of OPUS user interface


	Element	Definition
9	Status bar	The status bar informs about active tasks currently running in OPUS Introduction, and shows their status.
10	Instrument status	The instrument status is indicated by the OPUS status light () as follows: <ul style="list-style-type: none">• Gray: no spectrometer connected• Green: spectrometer connected• Yellow: warning• Red: error

Table 3.2: Elements of OPUS user interface

3.4 Access control

OPUS has an access control system which is shown in figure 3.2.




Figure 3.2: Access control in OPUS

3.4.1 User record

The user records in the OPUS user database comprise:

- User ID
- Password
- Operator name
- User group (Administrator, Operator)
- Signature categories
- Accessible workspaces

3.4.2 OPUS workspace

The workspaces comprise:

- Menu contents
- Toolbar contents
- Access rights within OPUS
- Other settings (e.g. file increment mode, display limits)

Detailed information on the user management e.g. to add, modify or remove users from the user database are described in the OPUS Reference Manual.

3.4.3 Data recording

In OPUS, details on the single manipulations and/or evaluations on spectra files are saved in the HISTORY data block (chapter 4.3).

In addition, the global audit trail logs all actions which have been performed in OPUS in connection with the user and signature management, as well as with the single spectrum file, within a defined period of time.

Next to the type of action the global audit trail also shows the user ID of the user who is responsible for this action. Details on the global audit trail are described in chapter 4.4.


4 Working with OPUS

4.1 Measuring a spectrum

Chapter 4.1 et seq. describes the individual steps to be performed or observed BEFORE, DURING or AFTER measuring spectra. Observe the order of the steps described.

The dialogs depend on the spectrometer configuration currently used. Therefore, the following dialogs may be slightly different from the ones displayed on your screen.

4.1.1 Connecting to spectrometer

1. On the *Measure* menu, select the *Optic Setup and Service* command, or click the  icon from the toolbar.
2. In the dialog that opens on the *Optical Bench* tab, check whether the spectrometer settings are OK. If yes, close the dialog. If no, continue with step 3.
3. From the *Configuration* drop-down list, select the particular spectrometer type.
4. Enter the spectrometer's IP address into the *Optical Bench URL* entry field.
5. Click the *Connect* button.
6. Wait until the connection to the spectrometer has been established.
 - *The spectrometer is connected as soon as the spectrometer name and firmware version are displayed (figure 4.3). If the connection has been successful, the OPUS status light is green.*







Figure 4.3: Spectrometer connecting data

7. Click the *Save Settings* button.

4.1.2 Setting up measurement parameters

1. On the *Measure* menu, select the *Measurement* command, or click the  icon from the toolbar.
2. On the dialog that opens, define the measurement parameters on the different tabs.
 - *To learn more about error diagnostics when setting measurement parameters refer to chapter 4.1.2.1. Details on the individual measurement parameters are described in the OPUS Reference Manual.*
3. Click the *Accept & Exit* button.

4.1.2.1 Error diagnostics when setting measurement parameters

Display in OPUS	Meaning
	The yellow exclamation mark indicates that the particular tab contains a warning for one or more measurement parameters.
	 Measurement still possible.

Example:

- Yellow marking in front of the drop-down list or in the entry field of a measurement parameter.

Result spectrum: Absorbance Save data from: 3100 cm-1 to: 1000 cm-1

- ☞ If you position the cursor on the drop-down list or entry field, a tooltip explains why a warning is available for the particular measurement parameter.

Result spectrum: Absorbance No background spectrum available! cm-1 to: 1000 cm-1 BMS,DTC, limits spectral range to > 4000 cm-1

Table 4.3: Error diagnostics when setting measurement parameters

Display in OPUS



Meaning

The warning symbols indicate that the particular tab contains one or more measurement parameters^a, which are missing or have not been correctly defined as they do not match the optic configuration.

i Measurement not possible.

Example:

- Red marking in front of the drop-down list or in the entry field of a measurement parameter.



- ☞ If you position the cursor on the drop-down list or entry field, a tooltip explains the missing or wrong measurement parameter.



Table 4.3: Error diagnostics when setting measurement parameters

- Some types of measurement parameters interact with other ones, i.e. if one measurement parameter is corrected, the warning indicated for other measurement parameters affected by the correction may no longer apply.

4.1.3 Calibrating interferogram peak

i The exact interferogram peak position must be determined and stored BEFORE the very first measurement. Once determined, the interferogram peak position must be determined only if hardware changes have arisen.

1. On the *Measure* menu, Select the *Advanced Measurement* command.
2. Click the *Check Signal* tab.
3. On the dialog that opens, activate the *Interferogram* option button.
4. Use the arrow buttons (figure 4.4) to move the scan region to the left or right if no interferogram peak is displayed.




Figure 4.4: Moving scan region

5. Click the *Save Peak Position* button to store the interferogram peak.
 - *The exact peak position is required to be able to perform Fourier Transformation.*
6. Click the *Accept & Exit* button.

4.1.4 Storing experiment file

1. On the *Measure* menu, select the *Advanced Measurement* command.
2. Click the *Advanced* tab.
3. On the dialog that opens, define the number of scans in the *Sample/Background Scan Time* entry fields (figure 4.5).





Figure 4.5: Defining number of scans

4. Define the path to automatically store the measuring data, e.g. *MEAS.
 - Use the  button to exactly select the storage path. If OPUS is to run in the validation mode, special requirements on the storage path must be fulfilled. Details are described in chapter 4.1.4.1.
5. Determine the data type for the result spectrum.
 - It is recommended to select 'Absorbance' as result spectrum data type. Figure 4.6 shows possible result spectrum data types.




Figure 4.6: Result spectrum data types

- *The data type of the result spectrum can be changed at any time.*
6. Click the Save button.
 7. On the dialog that opens, define a name for the experiment file and save this name.
 - *It is recommended to repeatedly use the experiment file. When working with macros, experiment files are essential.*

4.1.4.1 Storing path in the validation mode

When working in the validation mode, files have to be saved into a protected data pool. The term protected data pool¹ refers to public directories and sub-folders in the `C:\Users\Public\Public Documents\Bruker\OPUS_<version>` path.

Paths which do not exclusively use public directories and/or sub-folders are not accepted by OPUS, and the user will not be allowed to save spectrum files into these directories.

1. In the protected data pool files can be edited but not deleted nor overwritten. Programs other than OPUS have no access to the files stored in the protected data pool.

4.1.5 Measuring background spectrum

In case of a background spectrum you measure an open-beam spectrum without any sample in the optical path. Therefore, a background spectrum is also called single-channel reference spectrum. Before starting sample measurement you have to measure a background spectrum, if not yet available.

The background spectrum reflects the effects that the spectrometer itself has on measured spectra. For example, the performance of the light source varies at some frequencies, or the optics absorbs light stronger at a certain wavelength.

1. Make sure that the OPUS status light is green.
2. Make sure that there is no sample available in the sample compartment.
3. On the *Measure* menu, select the *Advanced Measurement* command.
4. Click the *Optic* tab.
5. On the dialog that opens, click the *Aperture* setting drop-down list and select the same value used to acquire a sample spectrum.
6. Click the *Basic* tab.
7. On the dialog that opens, click the *Background Single Channel* button.
 - *The status bar below the spectrum window indicates the progress of the background measurement, e.g. the number of scans already measured. Background measurement has finished if the status bar reads 'No Active Task'.*
8. Wait until the background measurement has finished, then measure the sample spectrum (chapter 4.1.6).

4.1.6 Measuring sample spectrum

1. If the background measurement (chapter 4.1.5) has finished, place the sample into the optical path of the spectrometer.
 - *How to do this depends on the spectrometer configuration.*
2. On the *Measure* menu, select the *Advanced Measurement* command.
3. Click the *Basic* tab.

4. On the dialog that opens, define the sample description and sample form in the particular entry field.
 - *This information is stored together with the spectrum.*
5. Click the *Sample Single Channel* button to start measurement.
 - *The 'Measurement' dialog closes and the spectrum window opens. The status bar below the spectrum window indicates the progress of the sample measurement. The sample measurement has finished if the measured spectrum is displayed in the spectrum window (figure 4.7).*




Figure 4.7: Measured sample spectrum displayed in the spectrum window

Strictly speaking, the spectrometer measures an interferogram rather than a spectrum. Using Fourier Transformation this interferogram is transferred into a single-channel spectrum. On the basis of this single-channel spectrum and the background spectrum measured, a result transmission spectrum is calculated.

To better understand the different types of spectra, figure 4.8 shows the single-channel spectrum, background spectrum and the calculated transmission spectrum. The particular data blocks are displayed in the browser.

Details on the data blocks are described in the OPUS Reference Manual.




Figure 4.8: All spectrum types of a measured sample spectrum

	Definition	Data block in the browser
1	Transmission spectrum	TR
2	Background spectrum (single-channel)	RSC
3	Sample spectrum (single-channel)	SSC

Table 4.4: Definition of spectrum types

4.2 Baseline correction

Directly after sample measurement has finished, the baseline correction is useful to manipulate spectra which could not have been measured very well due to the sample material available.

If the baseline of the spectrum is sloped, curved or significantly below 100% transmission, the sample preparation might have been insufficient. Instead of measuring the sample again, perform a baseline correction. In many cases, this will solve the problem.

However, you should always try to prepare the samples properly to get good spectra, instead of correcting them after measurement.

Figure 4.9 shows an example for baseline correction. The upper spectrum is the original spectrum, the lower one is baseline corrected.





Figure 4.9: Spectrum before (top) and after (bottom) baseline correction

4.2.1 Performing baseline correction

1. On the *Manipulate* menu, select the *Baseline correction* command, or click the  icon from the toolbar.
2. On the browser, double click the spectrum block (AB or TR) of the sample spectrum measured.
 - *The spectrum file is displayed in the 'File(s) to correct' selection field.*
3. Click the *Select Method* tab.
4. On the dialog that opens, select a baseline correction method (figure 4.10).
 - *Details on the single method types are described in the OPUS Reference Manual.*




Figure 4.10: Baseline correction - Select method

5. Enter the number of baseline points manually.
 - *It is recommended to accept the default value of 64.*
6. Activate the *Exclude CO₂ bands* checkbox to exclude the CO₂ range during baseline correction.
7. Click the *Correct* button.
 - *Baseline correction will be executed immediately.*

The browser window shows a red document symbol (figure 4.11) on top of the blue document symbol (original data). The red document symbol indicates that the spectrum file has been manipulated.




Figure 4.11: Spectrum file displayed in browser after baseline correction

4.2.1.1 Baseline correction completed

After baseline correction has been completed the spectrum has not yet been saved to hard disk and exists only as a temporary file. Save the file on hard disk, use a different file name to prevent the original file from being overwritten (chapter 4.5.2).

4.2.1.2 Undo baseline correction

1. On the browser, right click the red document symbol.
2. From the pop-up menu that opens select the *Undo all Manipulations* command.
3. On the dialog that opens, click the *Undo* button (figure 4.12).




Figure 4.12: Undo baseline correction

- *The red document symbol is no longer displayed in the browser, and the original data are available again.*

4.3 History of spectrum files


If a spectrum file loaded has been edited or changed or if global measuring settings have been changed, all these actions can be viewed in the non-editable, non-deletable history data block ().¹

Figure 4.13 shows parameters from a HISTOTRY data block of a spectrum file with several types of manipulations performed.

1. In addition, OPUS uses the global audit trail to report and log measurements, manipulations and/or evaluations of spectrum files, as well as changes made on particular setting options.

- ☞ To have the parameters displayed right click on the red document symbol and select the *Show Report* command.




Figure 4.13: Parameters of HISTORY data block

The HISTORY data block contains the following parameters:

- user ID
- OPUS version
- spectrum file name
- all previous manipulations (in the order of performance)

The parameters are stored in one single spectrum file, together with the spectrum edited. This ensures that all data manipulations are recorded, e.g. type of manipulation method, time etc.

4.4 Viewing global audit trail

i The global audit trail records all actions which have ever been taken in OPUS, in chronological order. Details on the global audit trail are described in the OPUS Reference Manual.

1. On the *Validation* menu, select the *View global audit trails* command.
2. In the view that opens, define a start and end date for the search period.
3. Click the *Request* button.
 - *The query result is displayed in chronological order in a table (figure 4.14).*

Activity	Object ID	Object snapshot	Time stamp	User ID	User full name
<all>	<all>	<all>	<all>	<all>	<all>
Login	32a60100-128...	Admin	2017/11/17 11:23:...	Admin	Admin
Logout	32a60100-128...	Admin	2017/11/17 12:46:...	Admin	Admin
Login	32a60100-128...	Admin	2017/11/17 12:46:...	Admin	Admin
Logout	32a60100-128...	Admin	2017/11/17 14:35:...	Admin	Admin
Login	32a60100-128...	Admin	2017/11/17 14:37:...	Admin	Admin
Logout	32a60100-128...	Admin	2017/11/17 14:40:...	Admin	Admin
Login	32a60100-128...	Admin	2017/11/20 07:54:...	Admin	Admin
Modify file		Abboe05.0	2017/11/20 12:13:...	Admin	Admin
Modify file		Abboe05.0	2017/11/20 12:14:...	Admin	Admin
Logout	32a60100-128...	Admin	2017/11/20 12:14:...	Admin	Admin
Login	32a60100-128...	Admin	2017/11/20 12:14:...	Admin	Admin
Logout	32a60100-128...	Admin	2017/11/20 14:54:...	Admin	Admin
Login	32a60100-128...	Admin	2017/11/21 07:00:...	Admin	Admin
Logout	32a60100-128...	Admin	2017/11/21 08:23:...	Admin	Admin
Login	32a60100-128...	Admin	2017/11/21 08:24:...	Admin	Admin

Figure 4.14: Example of a global audit trail result view

4.4.1 Loading HISTORY data block from global audit trail


i Modifications performed on a spectrum file are saved in the HISTORY data block, which is added to the spectrum file. The HISTORY data block can be directly loaded from the global audit trail for a particular spectrum file.

1. In the result view of the global audit trail, search the spectrum file which has been modified.
2. Double click the name of the modified spectrum file in the *Object Snapshot* column.
 - *The parameters of the HISTORY data block are displayed in a separate report view.*

4.5 Saving spectrum


i By default, the sample spectrum measured (original file) is saved in the path defined in chapter 4.1.4. If the original file has been manipulated, it is recommended to save the original file by using a different file name (chapter 4.5.2), to prevent the original file from being overwritten. Once overwritten, restoring the original file is not possible anymore.

4.5.1 Saving and overwriting the original file

1. On the *File* menu, select the *Save File* command, or click the  icon from the toolbar.
2. On the dialog that opens, drag & drop the spectrum file from the browser to the *File(s) to save* list box.

3. Click the *Save* button.
 - *The original file is overwritten. Restoring the original data is not possible anymore.*

4.5.2 Saving without overwriting the original file

1. On the *File* menu, select the *Save File As* command, or click the  icon from the toolbar.
2. On the dialog that opens, drag & drop the spectrum file from the browser to the *File to save (OPUS format)* list box (figure 4.15).





Figure 4.15: Save File As

3. Define the file name.
 - *By default, OPUS provides a file name which can always be changed.*
4. Define the path.
 - *The 'Change Path' button allows to exactly select the storage location.*

5. Activate the 'Overwrite' checkbox if a previous file is to be overwritten.
 - *If the checkbox is deactivated, the file name is incremented. The 'Overwrite' checkbox is disabled when the spectra have been measured in validation mode.*
6. Make the settings on the *Mode* and *Data Point Table/XML* tab.
 - *The spectrum file is normally saved in OPUS format. However, the 'Mode' and 'Data Point Table/XML' tabs allow to save the file in different formats, e.g. JCAMP-DX or XML format, as a plain XY data table, in Pirouette (DAT) or in GALACTIC format (SPC). Details on the tabs are described in the OPUS Reference Manual.*
7. Click the *Save* button.

4.6 Load spectrum file

1. On the *File* menu, click the *Load File* command, or click the  icon from the toolbar.
2. On the dialog that opens, select the particular spectrum file.
 - *The dialog is described in chapter 4.6.1.*
3. Click the *Open* button.
 - *The spectrum is displayed in the spectrum window.*

4.6.1 Description of *Load Spectrum* dialog




Figure 4.16: *Load Spectrum* dialog

	Element	Definition
1	Title bar	The title bar shows the directory in which the spectrum files are saved.
2	Spectrum files	Depending on the data type (4) selected all spectra files are displayed, which are available in the particular directory.
3	Preview	The preview shows the spectrum without axes, and the data blocks of the spectrum file.


Table 4.5: Elements of the *Load File* dialog

	Element	Definition
4	Data type	The drop-down list contains different data types. Spectra measured in OPUS have the <i>OPUS Spectrum</i> data type.
5	Parameter	By default, 3 parameters are displayed per each spectrum file. The drop-down lists contain different types of parameters. If a spectrum file has been selected, the values are displayed next to the parameters on the right (figure 4.16). As not all parameters must be specified for each spectrum, some of the parameter values may be missing even if the spectrum file has been selected.


Table 4.5: Elements of the *Load File* dialog

4.7 Printing spectrum file

4.7.1 Using the *Quick Print* command

- ☞ On the *Print* menu, select the *Quick Print* command, or click the  icon from the toolbar.
 - *Depending on the window type currently used, OPUS selects an appropriate plot layout. For example, if you work in the spectrum window, the spectrum will be printed within a frame. If you search a spectrum library, the results as well as a hit list, spectral information and the structure of the product will be printed. More details on pre-defined layout plots are described in the OPUS Reference Manual.*

4.7.2 Using the *Print Spectra* command

1. On the *Print* menu, select the *Print Spectra* command, or click the  icon from the toolbar.
2. On the dialog that opens (figure 4.17), you can individually define the basic settings for the print-out.




Figure 4.17: *Print Spectra* dialog

3. Select a layout.
4. If required, make further settings on the *Frequency Range* or *Options* tabs.
 - *More details on the single tabs are described in the OPUS Reference Manual.*
5. Click the *Print* button.


5 Analysis Report

OPUS allows to depict sample parameters and evaluation results for one or several spectrum files in an analysis report. The spectra files need to contain the data block of the respective evaluation. Evaluation data blocks are e.g. PEAK, INTEG, QUANT, IDENT, SEARCH, CT.



Figure 5.18: Spectrum file with data blocks of different evaluation types

5.1 Generating report

1. On the *Print* menu, select the *Generate Analysis Report* command, or click the  icon from the toolbar.
2. On the dialog that opens, click the *Load* button.
3. Select the particular spectrum file and click the *Open* button.
 - *The spectrum file is displayed in the 'File(s) for report generation' list box (figure 5.19).*




Figure 5.19: Generating analysis report

4. Define the report layout.
 - *Single sample report*¹: one report for each spectrum file
 - *Multi sample report*²: one report for all spectrum files
 5. Select the template file.
 - *By default, OPUS provides several types of template files. Template files have the extension *.art and are stored in the C:\Users\Public\Public Documents\Bruker\OPUS_<version> directory.*
 6. Click the *Output to* tab and define how the report is to be printed.
 - *If you want to print the report as PDF or csv file, you have to define the file name of the report and the path to store the report.*
 7. Click the *Generate* button.
-
1. To generate a single sample report you can use several spectra files at once. In this case, a separate report will be generated for each spectrum file separately.
 2. A multi sample report lists all spectrum files with the respective data in one report.

Index

A

Access control	17
Analysis report.....	43
Audit trail.....	18, 33
Loading HISTORY data block.....	36

B

Background spectrum	27, 29
Baseline correction	30, 31
Browser.....	15

C

Calibrating interferogram peak	24
--------------------------------------	----

E

Experiment file.....	25
----------------------	----

H

HISTORY data block	33, 34
--------------------------	--------

I

Installing	11
Instrument status	16

L

Login.....	12
------------	----

M

Measurement parameters	21
Menu bar	14
Multi sample report	44

O

Online help	15
OPUS view.....	15
Overview window	15

P

Preview	39
---------------	----

R

Result spectrum	26, 28
-----------------------	--------

S

Sample measurement	28
Sample spectrum	29
Service Hotline	6
Single sample report.....	44
Spectrum	
Loading	38
Printing	41
Saving	36
Spectrum window	15, 28
Status bar	16

T

Terms.....	5
Title bar	14
Toolbar	15

U

User ID.....	17
User interface	14

V

Validation mode.....	12, 13
----------------------	--------

